# Ambu® Neuroline Subdermal Needle Electrode

## Instructions for use

## Instrukcijas za polzvanе

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## Subdermal Needle Electrode

### Recommended Applications

Subdermal electrode has leadwire with touch proof safety connector DIN 42-802. Subdermal Electrode is designed for direct connection with IEC 60601-1 compliant Evoked Potentials (EP), Intra-Operative Monitoring (IOM) or other Neurophysiological equipment with DIN 42-802 input socket.

### Intended Use

The Ambu® Neuroline Subdermal Needle Electrode is intended for use with recording, monitoring and stimulation equipment for the purpose of recording of biopotential signals for example: Electromyography (EMG), Electroencephalograph (EEG) and Nerve potential signals. The electrodes are sterile and for single use only.

The Subdermal electrode should be used by health care professionals trained in performing neurophysiological examinations.

The Subdermal electrode is for use in hospital/clinic settings.

### Precautions

- For single use only. Reuse of the product can cause contamination, leading to infections.
- The needle electrode shall be administered to and removed from patients by manipulating the needle and not the leadwire.
- DO NOT soak, rinse, or sterilize this device as these procedures may leave harmful residues or cause malfunction of the device. The design and material used are not compatible with conventional cleaning and sterilization procedures.
- During surgical procedures, the electrode should be placed as far as possible from the equipment which generates electromagnetic interference (EMI) to minimize unwanted radiofrequency interference (RFI).
- Simultaneous use of needle electrode with a surgical equipment operating at high frequency may result in burns at the site where needle electrode is applied. The needle electrodes should be placed as far as possible from the electro-surgical area to minimise unwanted RF current flow.
- Do not use if individual pouch is damaged or has been opened.
- Do not reuse.
- Rx only.


### Sterilization/Sterilisation

Sterilized/Sterilised by Ethylene oxide sterilization/sterilisation (ETO).

### Packaging

One needle electrode in a pouch. Do not use if package is damaged.

### Directions for use

Ambu® is a registered trademark of Ambu A/S, Denmark. Ambu A/S is certified according to ISO 13485

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