

Ambu® Disposable Open Cuff Face Mask

DATASHEET

Ambu Disposable Open Cuff Face Masks are a range of anatomically shaped face masks that feature a soft, colour-coded cuff for quick size identification and establishment of a tight, comfortable seal.

The face masks also feature a transparent dome that allow for visualisation of the patient's condition.

The size 0 mask is designed for prematurely born infants and features a pliable, circular shape.

Available in 8 sizes, the range of Ambu Disposable Open Cuff Face Masks are suitable for use across infant, paediatric and adult patients.

These products are PVC-free and latex-free.



KEY BENEFITS

Lightweight anatomically shaped face masks with a soft cuff

Available in 8 different color-coded sizes for easy identification

Size 0 mask is designed for prematurely born infants and is circular in shape

Clear mask shell for visualisation of patient's condition

PVC-free and latex-free

Ambu

MATERIALS

Product	Material
Cuff	TPE
Dome	PolyPropylene
Premature size 0	Silicone

SIZE

Type	Weight
0 Premature	8.5 g
1 Neonate	7.1 g
2 Infant	8.9 g
3 Pediatric	15.1 g
4 Small Adult	18.5 g
5 Medium Adult	21.6 g
6 Large Adult	25 g
7 Extra Large Adult	27.8 g

ORDERING SPECIFICATIONS

Item No.	Type	Quantity per box
○ 306 0000 00	Ambu® Disposable Open Cuff Face Mask size 0 Premature	20
● 306 1000 00	Ambu® Disposable Open Cuff Face Mask size 1 Neonate	20
● 306 2000 00	Ambu® Disposable Open Cuff Face Mask size 2 Infant	20
○ 306 3000 00	Ambu® Disposable Open Cuff Face Mask size 3 Pediatric	50
● 306 4000 00	Ambu® Disposable Open Cuff Face Mask size 4 Small Adult	50
● 306 5000 00	Ambu® Disposable Open Cuff Face Mask size 5 Medium Adult	50
● 306 6000 00	Ambu® Disposable Open Cuff Face Mask size 6 Large Adult	50
● 306 7000 00	Ambu® Disposable Open Cuff Face Mask size 7 Extra Large Adult	20



Ambu® Disposable Open Cuff Face Mask



Ambu® Disposable Open Cuff Face Mask sizes 3, 4, 5, 6



Ambu® Disposable Open Cuff Face Masks sizes 0-7

Ambu

Distributed by:



 Baltorpbakken 13
 2750 Ballerup
 Denmark
 T +45 72 25 20 00
 ambu.com

Manufactured by:

Ningbo Shengyurui Medical Appliances Co., Ltd.
 No.138 Binhaisi Rd.,
 Hangzhou Bay New Zone
 315336 Ningbo City,
 Zhejiang Province, China


 0120

US: Rx Only