

# Ambu® Silicone Face Mask

## DATASHEET

The Ambu® Silicone Face Mask is designed for use with manual and automatic resuscitators, anaesthetic circuits and ventilators.

Ambu® Silicone Face Masks feature a flexible cushion manufactured from 100 % silicone that works to establish a tight yet comfortable seal. The cushion is self-inflating, and the pressure within the cushion may be regulated (except for size 0A). A transparent dome allows for a visual check of the patient's condition, and an integrated thumb rest supports hand positioning throughout ventilation.

Ambu® Silicone Face Masks can be reprocessed up to 30 times. They are suitable for both manual and automated cleaning and sterilization, and are capable of being autoclaved at a temperature of 134 °C (except size 0A). The flexible cushion is also removable to facilitate thorough cleaning.



## KEY BENEFITS

Suitable for reprocessing up to 30 times

MR Safe

Transparent dome for visual check of bleeding, vomitus and spontaneous breathing

Thumb rest to support hand positioning for a tight, easy seal

Comfortable design and texture

Flexible cushion manufactured from 100 % silicone which is also removable to allow for comprehensive cleaning

Suitable for autoclaving at a temperature of 134 °C/ 273.2 °F (except size 0A)

Available in 6 sizes

# Ambu

## ORDERING SPECIFICATIONS

Item no.	Type
000 313 000	Ambu® Silicone Face Mask size 0A
000 314 000	Ambu® Silicone Face Mask size 0 w.boring
000 315 000	Ambu® Silicone Face Mask size 0
000 316 000	Ambu® Silicone Face Mask size 2
000 312 000	Ambu® Silicone Face Mask size 3/4
000 317 000	Ambu® Silicone Face Mask size 5
000 319 000	Ambu® Silicone Face Mask size 6

## MATERIALS

Product	Material
Cuff	Silicone rubber
Cuff size 0A	Silicone rubber with foam filling
Dome	Polysulphone
Mask connector	Silicone rubber

## SIZE

Type	Weight
0A	28 g
0 w.boring	39 g
0	38 g
2	63 g
3/4	78 g
4	80 g
5	105 g
6	105 g



Ambu® Silicone Face Mask connected to Ambu® Mark IV Resuscitator.



Ambu® Silicone Face Mask connected to Ambu® Oval Plus Silicone resuscitators.

# Ambu

 **Ambu A/S**  
 Baltorpbakken 13  
 2750 Ballerup  
 Denmark  
 T +45 72 25 20 00  
 ambu.com

 <sub>2797</sub> US: Rx only