

## DATASHEET

The Intramuscular Injection Simulator is an innovative yet affordable simulator that is realistic in nature.

It fits on the area of the deltoid muscle as well as the vastus lateralis, rectus femoris, ventrogluteal, and dorsogluteal areas.

Because there is fluid in vessels within the simulated tissue, there is an actual possibility of seeing simulated blood in the hub of the needle on aspiration.

The skin and veins can be easily replaced, and the pad measures 17 cm x 8 cm x 4 cm.

# Intramuscular Injection Simulator



## KEY BENEFITS

- Innovative and affordable I.M. simulator
- Chance of aspiration
- Fits on deltoid muscle, vastus lateralis, rectus femoris, ventrogluteal and dorsogluteal areas
- Replaceable skin and veins

## ORDERING SPECIFICATIONS

### Item no.

LF01184U

### Intramuscular Injection Simulator, Light

incl. simulator, syringes, blood powder, hook-and-loop extension strap

# Ambu

📍 Ambu A/S  
Baltorpbakken 13  
2750 Ballerup  
Denmark  
T +45 72 25 20 00  
F +45 72 25 20 50  
ambu.com