

# Ambu<sup>®</sup> Man School



## Key features

- Ambu Hygienic system to prevent cross contamination
- Compact and portable training manikin
- Easy training of chest compressions
- Closed torso
- Easy placement of AED pads
- 10 manikins for easy class training

## Ambu Man School

Ambu Man School is a training manikin, which offers human appearance dedicated to the training of all procedures for basic life support resuscitation.

The manikin can be used for realistic training of external chest compressions and the placement of AED pads, in order to get used in the handling of AED's.

The unique hygienic system prevents cross contamination when demonstrating and training of the principals of ventilation through mouth and nose.

The Ambu Man School is intended to be used for education of pupils in schools under the supervision of an instructor/teacher.

### **Guideline conform**

The Ambu Man School allows to perform CPR training following the latest international CPR guidelines.

### **Compact and portable**

For an easy transportation & set up possibility the Ambu Man School has a low weight and compact packaging.

### **Easy AED training**

Due to the closed shaped torso, defibrillation pads can be placed, so simulated AED training can be performed realistically.

## Specifications

Weight	0,8 kg
Length	45 cm
Max. compression depth	6 cm

## Materials

Straps	Nylon
Base plate	Polycarbonate
Bellow	Thermoplastic Elastomer
Chest	Polypropylene
Head	Silicone
Face piece	Polyvinylchloride
Lung bag	Polyvinylchloride
Head bag	Polyethylene

## Order information

Item no.	Description
A268 007 000	Ambu Man School (Set of 10)
268 000 571	Face piece Ambu Man School (Set of 10)
268 000 572	Head bag Ambu Man School (Set of 100)
268 000 573	Accessories Ambu Man School (Set of 10)



Closed torso



Placement of AED pads



Ambu Man School hygienic system



Accessories

### Ambu A/S

Baltorpbakken 13  
DK - 2750 Ballerup  
Denmark  
Tel. +45 7225 2000  
Fax +45 7225 2053  
www.ambu.com